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Bilag:

Protocol for Impact genotyping – tagging, capture and amplification.

Kåre Lehmann Nielsen rev 02/2020

Input requirements: 50 ng highly purified gDNA A260/A280=1.8-2.0. (best option for quantification is qubit; however nanodrop required for ratio determination).

Abbreviations:

Qubit -Qubit 4 quick reference guide. ThermoFisher Scientific 2017.
NextProt -Nextera DNA library Prep reference guide (Document #15027987v01 – jan 2016) (Illumina)
MinElute -MinElute Reaction Cleanup kit using a microcentrifuge (Qiagen MinElute Handbook 01/2020)
HybProt -MyBaits hybridization capture for targeted NGS v4.01 April 2018
Tapestation -Agilent 2200 Tapestation system user manual 09/2015. Agilent technologies.

Per sample:

1. Follow NextProt tagmentation page 9
2. Following thermocycling end purify Rxns using MinElute.
3. Amplify the purified rxns using **a separate combination of barcodes** for each reaction for **8 cycles** according to the “Amplify Tagmented DNA “ section in NextProt page 12-14.
Note: Remember to note down which barcodes identifies which samples!
4. Pool 50% all reactions to a single tube
Note: store the remaining 50% in -20 freezer for backup.

Only a single pooled sample:

5. Purify Rxns using MinElute.
6. Quantify rxn using Qubit dsDNA HS assay, according to Qubit.
Note: Reference is not needed.
7. Analyze rxn using Tapestation HSD1000 following Tapestation page 56.
Note: use ladder if available.
8. Follow appendix A2 in HybProt page 15-16 using Phusion polymerase and P5reamp and P7 reamp primers.
9. Quantify rxn using Qubit dsDNA HS assay, according to Qubit.
Note: Reference is not needed.
10. Analyze rxn using Tapestation HSD1000 following Tapestation page 56.
11. Follow HybProt page 7-13 using Phusion polymerase – at page 12 split the rxn in three portions (5, 10, 10 μ l) and run 0, 8 and 14 cycles of PCR, respectively (3 reactions going into section 3.2).
12. Purify using MinElute.
13. Quantify rxn using Qubit dsDNA HS assay, according to Qubit.
Note: Reference is not needed.
14. Analyze rxn using Tapestation HSD1000 following Tapestation page 56.
Note: use ladder if available.
15. Determine A260/A260 by Nanodrop analysis.

Submit for sequencing using a single lane of HiSeq3000/4000 or similar throughput – about 300 mio reads are needed for an average of 250 reads pr site pr sample.